

🎯 ACOUSTIC WAVE

🎯 Cellulite treatment  
(PEFS stages I-III)

🎯 Improvement  
of skin elasticity

Body contouring

## WHAT IS THE ACOUSTIC WAVE

The treatment with AWS MED and SCULPTURE EVO devices is based on an innovative technology that uses acoustic vibrant waves to reduce the localized fat deposits and eliminate cellulite, with consequent improvement of the tone and elasticity of skin. The acoustic waves are propagated to the treated area releasing and transferring their energy to the fat cells.

The effect of this release of energy is the increase of bloodstream and lymphatic system. This process stimulates the synthesis of collagen, improving the hardening of the tissue.

The acoustic wave is generated by the gun shape probe, whose barrel is closed in a part by a metal cap against which is launched, using compressed air, a steel projectile, generating an acoustic wave (or sound wave) that is propagated from the point of contact with the skin to the deeper layers.


The high energy acoustic wave introduced into the body as high-pressure oscillations allow an effective and lasting treatment of cellulite (PEFS) stage I, II, III and not only.

## EFFECTS ON SKIN AND SUBCUTANEOUS TISSUE


Before treatment with acoustic waves: enlarged fat cells push upward towards the surface of the skin. The skin takes on the typical appearance orange peel. The inelastic vertical fibers of connective tissue do not stretch with the fat tissue around them. This gives to the skin the characteristic appearance as mattress. The local blood circulation is drastically reduced.

During the treatment with acoustic waves: high pressure acoustic waves are introduced into the body.

### Cellulite



Healthy skin



Skin with cellulite

⦿... This causes the relaxation of the connective tissue, and increases the tone of skin. The therapy improves the blood circulation into the tissue and stimulates the production of collagen, strengthening the skin.

⦿... After the treatment with acoustic waves: the skin and connective tissue are more toned and uniform finding their natural elasticity, with simultaneous visible and lasting improvement of the skin texture.



**E110142 - AWS MED COMPACT:**  
Desk device for acoustic wave treatments,  
certified as CE0476 medical device.



## APPLICATIONS IN AESTHETIC MEDICINE

- Body contouring
- Cellulite treatment (PEFS stages I-III)
- Firming of muscular and connective tissue
- Flattening of scar and wrinkles
- Improvement of skin elasticity
- Treatment of stretch marks due to pregnancy
- Circumference reduction

## AREAS TO BE TREATED

- Arms (triceps)
- Stomach
- Gluteus
- Hips
- Thighs



## ADVANTAGES

- Visible results after first sessions
- Technology is ideal for body-contouring and cellulite treatments
- Treatment suitable for any time of the year
- Long lasting effects
- No side effects



### **EI10141 -SCULPTURE EVO:**

Device for acoustic wave treatments,  
certified as CE aesthetic use device.

### **EI10140 - AWS MED:**

Device for acoustic wave treatments,  
certified as CE0476 medical device.



## Acoustic Wave

## Sculpture Evo AWS Med AWS Med compact

### Features and cutting-edge technology

- Color touch screen display
- Ergonomic handle, light and aesthetically enhanced
- One-Touch Mode
- 2 million shots guaranteed
- Pre-set protocols
- Replacing the actuator directly by the operator without the need for technical assistance
- Patient card to collect data, to develop a brief medical history and perform monitoring of treatments

### Technical sheet

	SCULPTURE EVO	AWS MED	AWS MED COMPACT
Main power supply	230 Vac, 50-60 Hz + 10%, (115 Vac, 50-60 Hz+10% - upon request)	230 Vac, 50-60 Hz + 10%, ((115 Vac, 50-60 Hz+10% - upon request)	230 Vac, 50-60 Hz + 10%, (115 Vac, 50-60 Hz+10% - upon request)
Display	10.4" touch screen color display	10.4" touch screen color display	6" touch screen color display
Shot emission modality	Burst auto, One touch, Burst, Continuous, Single	Burst auto, One touch, Burst, Continuous, Single	Continuous, Single, Burst
Work pressure	(1.5 - 5 ) bar 0.1 bar steps	(1.5 - 5 ) bar 0.1 bar steps	(1.5 - 4) bar 0.1 bar steps
Shot emission frequency	(1 - 20) Hz 1 Hz steps	(1 - 20) Hz 1 Hz steps	(1 - 15) Hz 1 Hz steps
Guaranteed shot number	2.000.000	2.000.000	2.000.000
Stored protocols	10 Woman -10 Men	Physiotherapy: 32 Aesthetic medicine: 11	Physiotherapy: 31 Aesthetic medicine: 11
Storable protocols on user's memory	200	200	200
Storable protocols on USB	4 Gb	4 Gb	Smart Card ( 250 protocols)
Dimensions	cm 62x39x111H	62x39x111H cm	39x28x17,5H cm
Weight	32 kg	32 kg	8 kg

### Accessori in dotazione

CODE	DESCRIPTION	SCULPTURE EVO	AWS MED	AWS MED COMPACT
ACC1317/2/E	USB key	●	●	●
ART1300/1	Pen for touch-screen display	●	●	●
ACC1336	Acoustic wave applicator with 35 mm transmitter	●	●	●
ACC1276	Interchangeable kit	●	●	
ART22315	35 mm transmitter			
ACC917	260 ml conductive gel			
ACC918	1000 ml conductive gel	●	●	●
ACC946	15 mm focused transmitter (AWS MED)		●	
ACC947	9 mm multi-focused transmitter (AWS MED)		●	
ACC948	15 mm multi-focused transmitter (AWS MED)		●	
ACC606/1	Kit of 10 Smart Cards (AWS MED Compact)			
ACC604	3-shelf trolley			
CONT72	Transport bag in TNT fabric			
	DVD Welcome Pack	●	●	●
	Smart Card			●